

Drug ordering instructions:

In EPIC: search TIMELESS in database. Complete the required fields. Add the note to pharmacy, highlighted below (choose applicable site location).

In the event that you are unable to use the EPIC order, follow these instructions to manually order study drug:

In the medication field, search "Non-formulary medication- INVESTIGATIONAL". Make sure that you are searching in the database and not facility list!

Complete the order with the following:

Drug name: "Tenecteplase 5mg/mL or Placebo Injection"

IRB #: "832114"

Study Name: "TIMELESS"

Dose: "___mg" (Fill in weight based down rounded to the nearest whole #)

Route: "Intravenous"

Frequency: "ONCE"

Admin Inst: "Administer via IV bolus over 5 seconds. MAX DOSE IS 25mg! Must be administered within 90 minutes of qualifying brain imaging."

Note to Pharmacy: "Drug and supplies are stored in the: HUP Central Inpatient Pharmacy/ PPMC Cupp central inpatient pharmacy/PAH IV Ante Room. Please refer to fax kit assignment for preparation. Dose should be administered ASAP after preparation, but reconstituted vials are stable for 8 hours, refrigerated."

Fax the IVRS assignment with kit number to inpatient pharmacy.

HUP: 215-349-5841

PAH: 215-829-5475

PPMC: 215-662-8218

Helpful tips for location of drug:

HUP pharmacy: Study drug is kept in the IDS cabinet, labeled TIMELESS P1831 and is non-refrigerated.

PAH pharmacy: Study drug is in a study bin in the IV Ante Room on the shelf next to the refrigerator.

PPMC pharmacy: Study drug is located in the Vaccine refrigerator in the central inpatient pharmacy, Cupp basement.